PENTAX Medical Portable Intubation Scopes

A Full Range of Easy-to-Use Portable Intubation Endoscopes





PENTAX Medical Portable Intubation Scopes A Full Range of Easy-to-Use Portable Intubation Endoscopes

For quick and reliable intubation and complete airway management, choose a PENTAX Medical compact, battery-powered, portable intubation endoscope. This modern system concept is ideal for use in the OR, ICU, or emergency room.

PENTAX MEDICAL QUALITY & SERVICE PROGRAMS

- The PENTAX Medical brand continues to be recognized for delivering cost-effective products of exceptional quality and durability.
- Industry-leading service plans, financing, technical service, and training.

PORTABILITY AND EASE-OF-USE

- Up to 160° tip deflection, a tight bending radius, and wide angle of view, simplifying orientation and the placement/positioning of endotracheal tubes.
- Ergonomically designed, lightweight control body, enhancing user performance.
- Battery-powered light source for true portability.
- PENTAX Medical CSV (Closed Suction Valve), a unique mechanism for improved suction.



• The perfect choice for OR, ICU, or emergency room situations.

SUPERB IMAGE QUALITY

- PENTAX Medical fiber technology delivers bright, clear, and high-resolution images.
- Upgradeable to video image through connection to an ocular-mounted camera system.

SPECIFICATIONS

Туре	Ultra-Slim	Slim	Standard	Large Channel	
Model Number	FI-7BS / FI-7RBS*	FI-9BS / FI-9RBS*	FI-10BS / FI-10RBS*	FI-13BS / FI-13RBS*	FI-16BS / FI-16RBS*
Angle of view (°)	95	90		95	
Field of view (°)			3-50		
Angulation (°) up/down	130/130		160/130		
Distal end (mm)	2.2	2.9	3.0	3.8	5.1
Insertion tube (mm)	2.4	3.1	3.5	4.2	5.2
Instrument channel (mm)		1.2	1.4	1.8	2.6
Working length (mm)	600				
Overall length (mm)	870	880			

*Available by special order

PENTAX CUSTOMER SERVICE: (800) 736-8293

PENTAX Medical 3 Paragon Drive Montvale, NJ 07645 Phone +1 800 431 5880 Fax +1 201 391 4189 www.pentaxmedical.com



Copyright 2014. All Rights Reserved. All company and product names and marks contained within are federally registered trademarks, trademarks or service marks of PENTAX of America, Inc. Product specifications are subject to change without notice, and without any obligation on the part of the manufacturer.