

# PENTAX Medical Fiber Bronchoscopes

## High-Resolution Fiber Bronchoscopy



**PENTAX**  
MEDICAL

# PENTAX Medical Fiber Bronchoscopes

## High-Resolution Fiber Bronchoscopy


*The highest optical quality in a complete product range.*

### PENTAX MEDICAL QUALITY & SERVICE PROGRAMS

- The PENTAX Medical brand continues to be recognized for delivering cost-effective products of exceptional quality and durability.
- Industry-leading service plans, financing, technical service, and training.

### FLEXIBILITY AND EASE-OF-USE

- Up to 180° tip deflection, a tight bending radius, and wide angle of view for easy maneuverability.
- Optimum size working channels (up to 3.2 mm) to accommodate a wide range of instruments and high suction capacity.
- PENTAX Medical CSV (Closed Suction Valve), a unique mechanism for improved suction.
- Compatible with a battery-powered light source (BS portable series only) or connected to a conventional light source.



- Ergonomically designed, lightweight control body, enhancing user performance.
- Fully immersible and reprocessible allowing these instruments to be completely and effectively high-level disinfected or sterilized.

### EXCELLENT IMAGE QUALITY

- PENTAX Medical fiber technology delivers bright, clear, high-resolution images for superior observation.

### SPECIFICATIONS

Type	Slim	Standard		Therapeutic		Portable	
Model Number	FB-8V	FB-10V	FB-15V	FB-18V	FB-19V	FB-15BS / FB-15RBS*	FB-18BS / FB-18RBS*
Angle of view (°)	100	120				100	
Focal range (mm)	2-50	3-50					
Angulation (°) up/down	180/130					180/130	180/90
Distal end (mm)	2.7	3.4	4.9	5.9	6.2	4.8	5.9
Insertion tube (mm)	2.8	3.5	4.9	6.0	6.2	4.9	6
Instrument channel (mm)	1.2	1.2	2.2	2.8	3.2	2.0	2.6
Working length (mm)	600						
Overall length (mm)	900					880	

\*Available by special order

**PENTAX CUSTOMER SERVICE: (800) 736-8293**

### PENTAX Medical

3 Paragon Drive  
Montvale, NJ 07645  
Phone +1 800 431 5880  
Fax +1 201 391 4189  
www.pentaxmedical.com

**PENTAX**  
**MEDICAL**