EB-1575K and EB-1975K Video Bronchoscopes



Next generation imaging for more efficient visualization

- New digital signal processing design to optimize colour and contrast
- Improved distribution of light to the peripheral area of the image for more effective use of field of view




Advanced imaging across all technology platforms

- Improved compatibility with PENTAX Medical HD+ platform
- Platform versatility choose the PENTAX
 Medical platform to meet your facility's needs

PENTAX Quality and Service Programs

- The PENTAX Medical brand continues to be recognized for delivering products of exceptional quality and durability
- Industry leading service plans, financing, technical service and training



Expanded reprocessing options

Sterrad[®] compatible



Easier intubation and access to target areas

- Full screen crisp, clear image
- Wide angulation range (210° for EB-1575K)
- Ergonomic control section for easy to reach suction valves and angulation lever

Lower cost per procedure

Reusable suction valve and inlet seal

EB-1575K / EB-1975K Specifications

Model Name	EB-1575K	EB-1975K
Direction of view	Forward	Forward
Field of view (°)	120	120
Tip Angulation up / down (°)	210/130	180/130
Ø Distal end (mm)	5.5	6.3
Ø Insertion tube (mm)	5.2	6.4
Ø Instrument channel (mm)	2.0	2.8
Working length (mm)	600	600
Total length (mm)	860	860

PENTAX Medical

1770 Argentia Road, Mississauga Ontario L5N 3S7 Canada Telephone Number: +1-905-286-5585 Fax Number: +1-905-286-5571 pentaxmedical.com

