


PENTAX Medical – Expand your possibilities for optical diagnosis.

The PENTAX Medical platform for high-end endoscopy.





Choose your system for success.

New possibilities for optical diagnosis with the unique premium series from PENTAX Medical.

High-Definition (HD+) endoscopy and optical diagnosis have become the new standard in clinical care. An increasing number of endoscopists regularly makes use of advanced technologies in their daily work. PENTAX Medical's imaging equipment can meet individual needs and technical requirements. The premium video processor EPK-i7000 and the high performance i10 series endoscopes offer the latest advances in endoscopic imaging technology.

PENTAX Medical provides the best solutions to support physicians all over the world in fulfilling their key roles in the health system. It provides aspiring endoscopists and leading physicians with highly sophisticated equipment and instruments to help sharing their latest findings with colleagues and fellows. Its innovative platform allows you to see better, achieve more and work more effectively.

See better

Exceptional image quality and full i-scan Management.


Achieve more

Double your information: more detailed information.

Work easier

With innovative features to support your daily challenges.





Angiodysplasia rectum




- HD+ high resolution imaging across the full range of endoscopes for better orientation
- Full i-scan Management to support detection, demarcation and complete pattern characterisation
- Close view of the mucosa for diagnostic as well as therapeutic procedures
- State-of-the-art image resolution (HD+) and i-scan functions combined in the Twin Mode for improved assessment of lesions
- Expanded therapeutic options in a standard endoscope: ratio outer diameter-working channel
- Videos, images and voice can be collected at the same time
- Integrated HD+ Video Recording and freeze scan to collect and share clinical findings
- Ergonomic design for improved comfort for users and patients
- Graduated flexibility (GDF) for a smooth intubation of the caecum




PENTAX Medical – The ultimate video processor: EPK-i7000.

Detect & define

Enhance your endoscopic options with HD+ and full i-scan Management.







Extraordinarily High-Definition (HD+) images and excellent illumination

- Supports fast orientation and detection
- Significant improvement in the visibility and evaluation of minute lesions
- Integrated zoom function for more detailed inspection of suspicious surface structures



- Detection and delineation support with i-scan Surface Enhancement (SE)
- i-scan SE retains the natural color tones
- Accentuation of tissue structures at the touch of a button
- Mucosal enhancement potentially supports the detection of flat lesions



Support in characterisation with i-scan Tone Enhancement (TE)


- Specific imaging technology for further assistance in endoscopic procedures
- Allows more accentuated display of mucosal structures which may support lesion characterisation
- Virtual chromoendoscopy may help to improve endoscopic diagnosis

Teach & assess

Twin Mode: State-of-the-art endoscopy with two images side-by-side.



The PENTAX Medical Twin Mode is an advanced technology that combines leading-edge HD+ quality and i-scan images to support detection and characterisation of lesions.



- Twin Mode doubles your information: different lesion details can be seen and compared at the same time
- Twin Mode allows simultaneous comparison of side-by-side endoscopic images during the entire examination
- Twin Mode facilitates the assessment of detected lesions by direct evaluation of different i-scan modes

Among tutors and speakers, the new Twin Mode is highly appreciated as an educational tool for "building the bridge" between HD+ white light images and i-scan findings.

The simultaneous comparison of enhanced clinical images is useful in teaching the appropriate characterisation of lesions.

Collect & share

Capture your examinations with integrated HD+ Video Recording.



It is all about collecting and sharing the latest scientific findings in the same HD+ image quality as experienced in the endoscopy room.





The EPK-i7000 provides integrated HD+ Video Recording and freeze scan to:

- Collect High-Definition (HD+) videos, images and voice
 Record your reports and personal comments on at the same time – at the touch of a button
- Share the most important findings within your scientific network
- the HD+ video
- Save cost: built-in Video Recording provides easy & quick recording

Diagnostic proof

Improved support along the diagnostic pathway.

For better detection

	HD+	HD+ and i-scan	p-value
n	931	1033	
Number of patients with an adenoma	27%	33%	P< .01
Number of adenomas per examined patient	.45	.59	P< .01
Number of non-adenomatous polyps/ examined patients	.51	.61	P=0.02
Number of advanced polyps/ examined patients	.061	.095	P< .01

The use of i-scan increases overall polyp detection rates, adenoma detection rates, and increases the number of advanced lesions detected during screening colonoscopy.

For better detection

Results	HD+	Standard*	p-value
Patients	100	100	
Patients with at least one adenoma	38	18	< .0001
Total number of lesions	145	37	< .0001
Total number of neoplasias	80	16	< .0001
Number of flat adenomas	22	3	< .0001
p<0.0001 Mann-Whitney Test			

With PENTAX HD+ endoscopes, the adenoma detection rate is significantly superior to standard resolution endoscopes.*

For easier characterisation

Pit pattern	Hyperplasia	Neoplasia	Total	Histology could be predicted with
Hyperplasia pit pattern I/II	63	2	65	a high sensitivity (98 %) and a high specificity (100 %) using i-scan TE.
Neoplastic lesions pit pattern III/IV/V	0	80	80	

[&]quot;High-definition colonoscopy combined with i-scan is superior in the detection of colorectal neoplasias compared to standard video colonoscopy." Hoffman A. et al, Endoscopy 2010

For more accurate diagnosis in GERD

	HD+	i-scan	Lugol's solution
Number of patients with macroscopically defined oesophagitis	9	12 ²	25¹
Number of small lesions	21	58¹	85
¹ Significant p<0.01 ² Not significant p<0.01			

i-scan as well as Lugol's solution (chromoendoscopy) significantly improve the identification of patients with oesophagitis.

For reliable surveillance in Barrett's oesophagus

i-scan	Histolog. Barrett	Histolog. Ø Barrett	Total
Type III-IV	26	2	28
Type I - II	2	13	15
Sensitivity: 92,8%, Specificity: 86,7%,			

Under the guidance of i-scan, targeted biopsies are reliable for predicting Barrett's epithelium.

A range that meets your individual needs.



EPK-i7000

The new premium HD+ video processor.

As fully equipped flagship model, the new EPK-i7000 video processor supports innovative and advanced endoscopic diagnosis. The EPK-i7000 is designed to provide individual solutions to master your daily challenges, to help with your research and to share your clinical findings.



i10 endoscope series

High-performance.

The PENTAX Medical i10 endoscope series delivers images of extraordinarily high resolution thanks to a high performance CCD chip. It enables you to fully exploit HD+ performance, offers expanded therapeutic options and gives a close view of mucosal details. The i10 series is the natural choice for endoscopic procedures.



EPK-i5000

From routine applications to challenging examinations.

The PENTAX Medical EPK-i5000 combines excellent illumination and HD+ quality at an optimal price-performance ratio. Ideal for practitioners who require routine functions offered by proven imaging technology in their day-to-day clinical work.



90i and 90K endoscope series

The redefinition of High-Definition.

The PENTAX Medical 90i and 90K endoscope series are the optimal choice for an exceptional high resolution image. Combined with excellent illumination, they support better orientation and detection in endoscopic procedures. With the 90i endoscope series you can fully experience HD+ during routine examinations.

This brochure is intended for distribution in EMEA.

[&]quot;High-Definition Colonoscopy Combined With i-scan Imaging Technology is Superior in the Detection of Adenomas and Advanced Lesions Compared to High-Definition Colonoscopy Alone."
Gagovic V. et al. Abstract DDW 2011

[&]quot;High-Definition Colonoscopy Combined With i-scan Imaging Technology is Superior in the Detection of Adenomas and Advanced Lesions Compared to High-Definition Colonoscopy Alone." Gagovic V. et al, Abstract DDW, 2011

[&]quot;High-definition endoscopy with i-scan and Lugol's solution for more precise detection of mucosal breaks in patients with reflux symptoms." Hoffman A. et al, Endoscopy 2009

[&]quot;Acetic acid guided biopsies compared with high-definition endoscopy and i-scan guided biopsies in the detection of Barrett's oesophagus." Hoffman A. et al. Poster Presentation, DDW 2010

PENTAX Europe GmbH

Julius-Vosseler-Straße 104 22527 Hamburg Germany

Tel.: +49 40 / 5 61 92 - 0 Fax: +49 40 / 5 60 42 13 E-mail: medical@pentax.de

PENTAX U.K. Limited

PENTAX House Heron Drive, Langley Slough SL3 8PN United Kingdom

Tel.: +44 17 53 / 79 27 33 Fax: +44 17 53 / 79 27 94 E-mail: medical@pentax.co.uk

PENTAX France Life Care S.A.S.

112 quai de Bezons B.P. 204 95106 ARGENTEUIL CEDEX France

Tel.: +33 1 / 30 25 94 78 Fax: +33 1 / 30 25 74 45 E-mail: contact.medical@pentax.fr

PENTAX Nederland B.V.

Amsterdamseweg 29 1420 CA Uithoorn Netherlands

Tel.: +31 88 / 5 30 30 30 Fax: +31 88 / 5 30 30 40 E-mail: medical@pentax.nl

PENTAX Italia S.r.l.

Via Dione Cassio, 15 20138 Milano

Italy

Tel.: +39 / 02 50 99 58 1 Fax: +39 / 02 50 99 58 60

E-mail: marketing.lifecare@pentaxitalia.it

SIMMEDICA – Sistemas Integrales de Medicina, S.A.

Avenida del Sistema Solar 25 28830 San Fernando de Henares · Madrid Spain

Tel.: +34 91 / 301 62 40 Fax: +34 91 / 751 31 15 E-mail: sim@simmedica.com

PENTAX Europe GmbH

Representative office in Moscow Sadovnicheskaya str, 82, build 2, entrance 6 Regus, Business center "Aurora"

Office 2012, 2013 Moscow, 115035 Russian Federation Tel.: +7 495 / 792 52 00 Fax: +7 495 / 792 35 66

PENTAX Turkey

Veko Giz Plaza, Meydan Sokak No:3/43 343396 Maslak – Istanbul

Turkey

Tel.: +90 212 / 705 05 26 Fax: +90 212 / 705 05 00

HOYA Corporation PENTAX Life Care Division

1-1-110, Tsutsujigaoka Akishima-shi 196-0012 Tokyo Japan

Tel.: +81 3 / 39 60 51 55 Fax: +81 3 / 53 92 67 24

