BROADEN YOUR PERSPECTIVE

Ultrasound Video Gastroscope EUS-J10 line-up



The PENTAX Medical Ultrasound Video Gastroscope EUS-J10 line-up is designed for optimized procedural efficiency. PENTAX Medical provides brilliant imaging and unique features in one clinical solution. This combination allows for clear visualization, successfully enabling endosonographic diagnosis and advancing therapeutic procedures for the benefit of patients.



Optimized procedural efficiency



Optimize access

Increased maneuverability and improved ergonomics contribute superb stability and ease of intubation to support higher procedural success rates, optimized working comfort and better patient tolerance. *The new way of doing.*

Enrich diagnostic findings

Greater detail, better delineation, and superior penetration: a proven PENTAX Medical image quality enhanced to provide improved diagnostic outcomes, and guidance on therapy. The new way of seeing



Redefine therapy

A unique, enlarged 4.0 mm working channel expands therapeutic capabilities to advance EUS-guided intervention while providing full control of the procedure. *The new way of treating patients.*

#InspiredbyPENTAXMedicalTripleAim

Ultrasound Video Gastroscope EUS-J10 line-up

Direction of view [°]	EG36-J10UR Forward	EG34-J10U Forward oblique (45°)	EG38-J10UT Forward oblique (45°)
Depth of field [mm]	3-100	3–100	3–100
Tip angulation up – down [°]	150/70	160/130	160/130
Tip deflection right – left [°]	70/70	120/120	120/120
Distal end width [mm]	10.4	12.9	14.3
Insertion tube width [mm]	12.1	11.6	12.8
Minimum instrument channel width [mm]	2.4	2.8	4.0
Insertion tube working length [mm]	1,250	1,250	1,250
Total length [mm]	1,566	1,566	1,566
Frequency [MHz]	5–13	5 –13	5 –13
Scan method	Radial	Convex	Convex
Scan angle [°]	360	150	150

EC REP G

EMEA Headquarter Germany

PENTAX Europe GmbH
Julius-Vosseler-Straße 104, 22527 Hamburg
T +49 (0)40 561 92-0 • F +49 (0)40 561 92-247
info.emea@pentaxmedical.com • www.pentaxmedical.com



6-10-1 Nishi-shinjuku Shinjuku-ku, Tokyo 160-0023 Japan

