


Gastrointestinal Endoscopes

Delivering innovation and performance for quality patient care



Gastrointestinal Endoscopes



Model Number	Product Description	Channel Diameter (mm)	Outer Diameter (mm)	Working Length (mm)	Up/Down Angulation (°)	Left/Right Angulation (°)	Angle of View (°)	Depth of Field (mm)	Forward Water Jet	Image Type	CCD Type	Control Body
EGD												
EG34-i10	Large Channel HD+	3.8	11.6	1050	210/120	120/120	140	2-100 (CloseFocus™)	Y	Video	HD+	i10 Ergonomic Control Body
EG29-i10	Standard HD+	3.2	9.8	1050	210/120	120/120	140	2-100 (CloseFocus™)	Y	Video	HD+	i10 Ergonomic Control Body
EG27-i10	Slim HD+	2.8	9.0	1050	210/120	120/120	140	2-100 (CloseFocus™)	-	Video	HD+	i10 Ergonomic Control Body
EG16-K10	Small Caliber	2.0	5.4	1100	210/120	120/120	140	4-100	-	Video	SD	i10 Ergonomic Control Body
EG-3890TK	Dual Channel	3.8 & 2.8	12.8	1050	180/120	120/120	140	4-100	Y	Video	SD	90 Series
EG-2990K	Standard SD	2.8	9.8	1050	210/120	120/120	140	4-100	Y	Video	SD	90 Series
EG-2790K	Slim SD	2.8	9.0	1050	210/120	120/120	140	4-100	-	Video	SD	90 Series
EG-2490K	8mm Slim	2.4	8.0	1050	210/120	120/120	140	4-100	-	Video	SD	90 Series
ERCP												
ED34-i10T2	DEC™ HD Duodenoscope	4.2	11.6	1250	120/90	90/105	100 (10 retro)	4-60	-	Video	HD	i10 Ergonomic Control Body
ED34-i10T	HD Duodenoscope	4.2	11.6	1250	120/90	90/105	100 (10 retro)	4-60	-	Video	HD	i10 Ergonomic Control Body
EUS (Endoscopic Ultrasound)												
EG-3870UTK	Linear Ultrasound Endoscope	3.8	12.8	1250	130/130	120/120	120	5-100	-	Video/US	SD	70 Series
EG-3270UK	Slim Linear Ultrasound Endoscope	2.8	10.8	1250	130/130	120/120	120	5-100	-	Video/US	SD	70 Series
EG-3670URK	Radial Ultrasound Endoscope	2.4	12.1	1250	130/60	60/60	140	4-100	-	Video/US	SD	70 Series
Small Bowel Enteroscopy												
VSB-2990i	Small Bowel HD Enteroscope	2.8	9.8	2200	180/180	160/160	140	3-100	-	Video	HD	90 Series
Colonoscopy												
EC38-i10L	Standard HD+	3.8	13.2	1700	180/180	160/160	140	4-100	Y	Video	HD+	i10 Ergonomic Control Body
EC38-i10NL	Standard HD+ w/ SCOPEPILOT®	3.8	13.2	1700	180/180	160/160	140	2-100 (CloseFocus™)	Y	Video/Nav	HD+	i10 Ergonomic Control Body
EC34-i10L	Slim HD+	3.8	11.6	1700	180/180	160/160	140	2-100 (CloseFocus™)	Y	Video	HD+	i10 Ergonomic Control Body
EC34-i10TL	RetroView™ HD+	3.2	11.6	1700	210/180	160/160	140	2-100 (CloseFocus™)	Y	Video	HD+	i10 Ergonomic Control Body
EC-2990Li	UltraSlim HD	2.8	9.8	1700	180/180	160/160	140	3-100	Y	Video	HD	90 Series
EC-3890TLK	Dual Channel	3.8 & 2.8	13.2	1700	180/180	160/160	140	3-100	Y	Video	SD	90 Series
EC-3890LK	Standard SD	4.2	13.2	1700	180/180	160/160	140	3-100	Y	Video	SD	90 Series
EC-3490LK	Slim SD	3.8	11.6	1700	180/180	160/160	140	3-100	Y	Video	SD	90 Series

PENTAX Medical

6715 Millcreek Drive, Unit 1
Mississauga, ON L5N 5V2
Phone +1 800 750 5558
Fax +1 905 286 5571
www.pentaxmedical.ca

PENTAX
MEDICAL
Excellence in Focus