




PENTAX Medical endoPRO iQ[®] Technology Solutions

A Foundation of Efficient Patient Care And Streamlined Operations



A comprehensive information and image management system that supports all stages of your endoscopy lab's clinical workflow, from patient arrival to discharge



PENTAX
MEDICAL

endoPRO iQ® Modules At A Glance

Functionality	Details	endoPRO iQ® Base Platform	endoPRO iQ® Video Recording/Editing Platform	endoPRO iQ® Advanced Documentation Platform	endoPRO iQ® Advanced Integration Platform
Scheduler	<ul style="list-style-type: none"> Unique whiteboard scheduler to view lab status at a glance Intuitive click, drag and drop design Customizable, color coded feature to easily track patient statuses Easy navigation to any date in the past, present or future Multiple facility/department support Optional worklist view 	X			
Patient/Procedure Data Entry	<ul style="list-style-type: none"> Easy one click access to all data entry screens Comprehensive patient data screen to easily enter patient demographics, allergies, medical history, current medications, family history and insurance Easy data entry for relevant procedure data such as indications, sedation plan, planned scopes, planned procedures, assistants, attending fellow and referring MD Reliable data accuracy; all data captured flows to downstream physician and nursing reports Access to all previous procedures, along with procedure notes and other scanned documents 	X			
Patient/Procedure Details Recall	<ul style="list-style-type: none"> Quick recall of all patient and procedure data by patient name, MR or SSN with instant access to all images, reports and videos Easy review of patient history with one click 	X			
Database Reporting	<ul style="list-style-type: none"> Instant access to comprehensive list of pre-defined database reports Instant access to custom data report generator for research projects and compliance requirements Automated collection of lab efficiency and equipment utilization data 	X			
Still Image Capture and Management	<ul style="list-style-type: none"> Digital image capture and storage from scope buttons in SD or HD (on endoscopic processors) Instant recall of captured images for review Image annotation tools to mark images with location and finding Extensive image editing and enhancement tools such as exposure and color correction Ability to print image reports for patient consultations 	X			
Video Capture and Storage	<ul style="list-style-type: none"> Digital video capture and storage from scope buttons in SD or HD (on endoscopic processors) Ability to capture unlimited video duration and frequency Instant archiving of digital video with the patient record Instant recall video clips for review 		X		
Video Editing	<ul style="list-style-type: none"> Instant and intuitive video editing to create an "executive summary" Still image capture from live video clips Video clip indexing with captions describing what is being shown Video clip tagging with searchable terms 		X		
Video Export	<ul style="list-style-type: none"> Video clip exporting to standard Windows supported format Video clip sharing and collaboration with surgeons when referring patients Video clip exporting for teaching or presentation purposes 		X		
Doc-U-Scribe™ Advanced Procedure Note Writer	<ul style="list-style-type: none"> Clinically relevant medical content developed in conjunction with leading academic institutions Support for most common endoscopy procedures, including double procedures Quick and intuitive operative note generation in minimal time Sophisticated language engine to create more natural-sounding reports Embedded color images in operative reports Participation in GIQuIC outcomes registry 			X	
CPT/ICD Support	<ul style="list-style-type: none"> Summary of suggested CPT codes on operative notes Summary of suggested diagnosis and history codes on operative note Tools to assist with faster reimbursements 			X	
Nursing Notes Advanced Nursing Documentation System	<ul style="list-style-type: none"> SGNA compliant for nursing documentation Quick charting for all phases of the patient encounter Easy vitals and medication totals tracking Quick and simple pathology module for tracking samples obtained, printing labels and lab requisitions 			X	
H7 Inbound Interface	<ul style="list-style-type: none"> Interface with hospital scheduling system to receive ADT/SIU/ORM messages Automated incorporation schedule data directly into the endoPRO system 				X
H7 Outbound Interface (requires Doc-U-Scribe™ Reporting Platform)	<ul style="list-style-type: none"> Outbound Results (ORU) interface to facility HL7 based system Seamless submission of selected data from operative notes to all major HIS systems such as MEDITECH, Cerner, EPIC, McKesson Star and Veterans Administration VistA and Clinical Procedures 				X
HL7 PDF Option (requires Doc-U-Scribe™ Reporting Platform)	<ul style="list-style-type: none"> Automatic export of link to formatted procedure report in PDF format for EMRs that support this functionality PDF copies contain embedded endoscopic images 				X
Vitals Monitor Integration (requires Nursing Notes Reporting Platform)	<ul style="list-style-type: none"> Automated capture from supported vitals monitors to seamlessly integrate vitals data into nursing report Integration with selected Data Scope, GE and Phillips Monitors 				X

endoPRO iQ® 7.6 is built on the latest technologies including Windows 7 and Windows Server 2008 client-server technologies powered by Microsoft SQL Server relational database management system.

PENTAX Medical

3 Paragon Drive
 Montvale, NJ 07645
 Phone +1 800 431 5880
 Fax +1 201 391 4189
www.pentaxmedical.com

