


# Gastrointestinal Endoscopes

*Delivering innovation and performance for quality patient care*



# Gastrointestinal Endoscopes



Model Number	Product Description	Channel Diameter (mm)	Insertion Tube Diameter (mm)	Working Length (mm)	Up/Down Angulation (°)	Left/Right Angulation (°)	Field of View (°)	Depth of Field (mm)	Forward Water Jet	Image Type	Sensor Type	Control Body
<b>EGD</b>												
OPTIVISTA plus Endoscopy Platform												
EG34-i10	Large Channel HD+	3.8	11.6	1050	210/120	120/120	140	2-100 (CloseFocus™)	Y	Video	HD+ CCD	i10 Ergonomic Control Body
EG29-i10	Standard HD+	3.2	9.8	1050	210/120	120/120	140	2-100 (CloseFocus™)	Y	Video	HD+ CCD	i10 Ergonomic Control Body
EG27-i10	Slim HD+	2.8	9.0	1050	210/120	120/120	140	2-100 (CloseFocus™)	-	Video	HD+ CCD	i10 Ergonomic Control Body
EG-2990i	Standard HD+	2.8	9.8	1050	210/120	120/120	140	5-100	Y	Video	HD+ CCD	90 Series
EG-2990Zi	MagniView HD	2.8	9.8	1050	210/120	120/120	140	Normal (WIDE) 5-100	Y	Video	HD CCD	90 Series
EG16-K10	Small Caliber SD	2.0	5.4	1100	210/120	120/120	140	4-100	-	Video	SD CCD	i10 Ergonomic Control Body
EG-1690K	Transnasal SD	2.0	5.4	1100	210/120	120/120	120	4-100	-	Video	SD CCD	90 Series
EG-3490K	Large Channel SD	3.8	11.6	1050	210/120	120/120	140	4-100	Y	Video	SD CCD	90 Series
EG-3890TK	Dual Channel SD	3.8 & 2.8	12.8	1050	180/120	120/120	140	4-100	Y	Video	SD CCD	90 Series
EG-2990K	Standard SD	2.8	9.8	1050	210/120	120/120	140	4-100	Y	Video	SD CCD	90 Series
EG-2790K	Slim SD	2.8	9.0	1050	210/120	120/120	140	4-100	-	Video	SD CCD	90 Series
EG-2490K	8mm Slim SD	2.4	8.0	1050	210/120	120/120	140	4-100	-	Video	SD CCD	90 Series
EG-290Kp	Standard SD	2.8	9.8	1050	210/120	120/120	140	5-100	-	Video	SD CCD	90 Series
IMAGINA Endoscopy Platform												
EG29-i10c	Standard HD+	3.2	9.8	1050	210/120	120/120	140	3-100	Y	Video	HD+ CMOS	i10 Ergonomic Control Body
VERSA Endoscopy Platform												
EG29-V10c	Standard HD	2.8	9.6	1050	210/90	100/100	140	3-100	-	Video	HD CMOS	V10 Series
EG27-V10c	Standard HD	2.8	9.0	1050	210/90	100/100	140	3-100	-	Video	HD CMOS	V10 Series
<b>ERCP</b>												
OPTIVISTA plus Endoscopy Platform												
ED34-i10T2	DEC™ HD Duodenoscope	4.2	11.6	1250	120/90	90/105	100 (10 retro)	4-60	-	Video	HD CCD	i10 Ergonomic Control Body
ED34-i10T	HD Duodenoscope	4.2	11.6	1250	120/90	90/105	100 (10 retro)	4-60	-	Video	HD CCD	i10 Ergonomic Control Body
ED-3490TK	SD Duodenoscope	4.2	11.6	1250	120/90	90/105	100 (10 retro)	4-60	-	Video	SD CCD	90 Series

Model Number	Product Description	Channel Diameter (mm)	Insertion Tube Diameter (mm)	Working Length (mm)	Up/Down Angulation (°)	Left/Right Angulation (°)	Field of View (°)	Depth of Field (mm)	Forward Water Jet	Image Type	Sensor Type	Control Body
<b>EUS (Endoscopic Ultrasound)</b>												
OPTIVISTA plus Endoscopy Platform												
EG-3870UTK	Linear Ultrasound Endoscope	3.8	12.8	1250	130/130	120/120	120	5–100	-	Video/US	SD	70 Series
EG-3270UK	Slim Linear Ultrasound Endoscope	2.8	10.8	1250	130/130	120/120	120	5–100	-	Video/US	SD	70 Series
EG-3670URK	Radial Ultrasound Endoscope	2.4	12.1	1250	130/60	60/60	140	4–100	-	Video/US	SD	70 Series
<b>Small Bowel Enteroscopy</b>												
OPTIVISTA plus Endoscopy Platform												
VSB-2990i	Small Bowel Endoscope HD	2.8	9.8	2200	180/180	160/160	140	3–100	-	Video	HD	90 Series
VSB-3430K	Small Bowel Endoscope SD	3.8	11.6	2200	180/180	160/160	140	3–100	-	Video	SD	30 Series
<b>Colonoscopy</b>												
OPTIVISTA plus Endoscopy Platform												
EC38-i10L	Standard HD+	3.8	13.2	1700	180/180	160/160	140	4–100	Y	Video	HD+ CCD	i10 Ergonomic Control Body
EC38-i10NL	Standard HD+	3.8	13.2	1700	180/180	160/160	140	2–100 (CloseFocus™)	Y	Video	HD+ CCD	i10 Ergonomic Control Body
EC-3890Li	Standard HD+	3.8	13.2	1700	180/180	160/160	140	4–100	Y	Video	HD+ CCD	90i Series
EC-3890LZi	MagniView HD	3.8	13.2	1700	180/180	160/160	140	Normal (WIDE) 4–100	Y	Video	HD CCD	90i Series
EC34-i10L	Slim HD+	3.8	11.6	1700	180/180	160/160	140	2–100 (CloseFocus™)	Y	Video	HD+ CCD	i10 Ergonomic Control Body
EC34-i10TL	RetroView™ HD+	3.2	11.6	1700	210/180	160/160	140	2–100 (CloseFocus™)	Y	Video	HD+ CCD	i10 Ergonomic Control Body
EC-3490TLi	RetroView™ HD	3.2	11.6	1700	210/180	160/160	140	3–100	Y	Video	HD CCD	90i Series
EC-2990Li	UltraSlim HD	2.8	9.8	1700	180/180	160/160	140	3–100	Y	Video	HD CCD	90 Series
EC-3890TLK	Dual Channel SD	3.8 & 2.8	13.2	1700	180/180	160/160	140	3–100	Y	Video	SD CCD	90 Series
EC-3890LK	Standard SD	4.2	13.2	1700	180/180	160/160	140	3–100	Y	Video	SD CCD	90 Series
EC-3490LK	Slim SD	3.8	11.6	1700	180/180	160/160	140	3–100	Y	Video	SD CCD	90 Series
EC-380LKp	Standard SD	3.8	13.2	1700	180/180	160/160	140	3–100	-	Video	SD CCD	80 Series
IMAGINA Endoscopy Platform												
EC38-i10cL	Standard HD+	3.8	13.2	1700	180/180	160/160	140	3–100	Y	Video	HD+ CMOS	i10 Ergonomic Control Body
EC34-i10cL	Slim HD+	3.8	11.6	1700	180/180	160/160	140	3–100	Y	Video	HD+ CCD	i10 Ergonomic Control Body
VERSA Endoscopy Platform												
EC38-V10cL	Standard HD	3.7	12.8	1650	180/180	160/160	140	3–100	-	Video	HD CMOS	V10 Series
EC38-V10cM	Standard HD	3.7	12.8	1300	180/180	160/160	140	3–100	-	Video	HD CMOS	V10 Series
<b>Sigmoidoscopy</b>												
OPTIVISTA plus Endoscopy Platform												
ES-3870K	Sigmoidoscope SD	4.2	12.8	700	180/180	160/160	140	3-100	Y	Video	SD	70 Series

**PENTAX Medical**

Americas Head Office  
3 Paragon Drive  
Montvale, NJ 07645  
Phone +1 800 431 5880  
Fax +1 201 391 4189  
[www.pentaxmedical.com](http://www.pentaxmedical.com)

