

PENTAX



ECY-1575K Flexible Video Cystoscope
PENTAX HD Platform and STERRAD® Compatibility

The New Standard in Flexible Video Cystoscopy and Reprocessing



ECY-1575K Specifications:

Model Name	ECY-1575K
Field of View (°)	120
Focal Range (mm)	3-50
Tip Angulation up / down (°)	210/130
Ø Distal End (mm)	5.5
Ø Insertion Tube (mm)	5.2
Ø Instrument Channel (mm)	2.0
Working Length (mm)	400
Overall Length (mm)	660

PENTAX Quality and Service Programs:

- > The PENTAX brand continues to be recognized for delivering products of exceptional quality and durability.
- > Industry leading service plans, financing, technical service and training.

Next Generation Imaging for More Efficient Visualization

- > New digital signal processing design to optimize colour and contrast.
- > Improved distribution of light to the peripheral area for more effective use of the field of view.
- > Enhanced clarity and brightness on a full-screen image.

User-friendly Design

- > Flexible angulation and tighter bending radius enable easier observation of the bladder and bladder neck.
- > Ergonomically-designed control body provides optimum physician comfort.


Advanced Imaging Across All Technology Platforms

- > Improved compatibility with PENTAX HD platform.
- > Platform versatility - choose the PENTAX platform to meet your facility's needs.

Expanded Reprocessing Options

- > STERRAD® compatible.

Objective Lens



Dual light guides

Instrument Channel



PENTAX Customer Service: (800) 750-5558