


## Single Use Therapeutic Accessories

*The Choice for Quality and Efficiency in Endoscopy.  
Available to ensure highest standards.*



**PENTAX Europe GmbH**  
Julius Vosseler Strasse 104  
22527 Hamburg  
Germany  
Tel: +49-40-561-920  
Fax: +49-40-560-4213

**PENTAX Medical Company**  
A Division of PENTAX of America, inc.  
3 Paragon Drive Montvale  
New Jersey 07645-1782  
USA  
Tel: +1-201-571-2300

**PENTAX Medical (Shanghai) Co., Ltd.**  
Room 701  
291 Fumin Road, Shanghai  
China 200031  
Tel: +86-21-6170-1555  
FAX: +86-21-6170-1655

**HOYA Corporation**  
PENTAX Lifecare Division  
1-1-110 Tsutsujigaoka  
Akishima-shi  
Tokyo 196-0012  
Japan



Product availability cannot be guaranteed in all countries. For further information please contact your local sales representative.



# PENTAX Medical

State-of-the-art accessories product line for diagnosis and treatment in endoscopy.



## Our commitment to quality

Safe and successful endoscopy depends upon the use of a finely tuned endoscopic system of components. This is the underlying principle behind PENTAX Medical's commitment to a product line of therapeutic accessories with the highest possible level of quality. We invite you to examine our carefully selected and meticulously crafted therapeutic accessories for state-of-the-art diagnostic and treatment applications. All PENTAX Medical accessories have been tested in conjunction with PENTAX Medical endoscopes for function and compatibility.

## How to select the correct accessories.

Before you select new accessories, please consider the following:


- The compatibility of your endoscope with the desired diagnostic/treatment application.
- The diameter of the instrument channel and working length of your endoscope.
- Use only accessories that are compatible with your endoscope to ensure the highest standards of endoscopic therapy.

## Single Use Biopsy Forceps




Description	Model Name	Compatible Endoscopes	
		Instrument Channel Diameter [mm]	Working Length [mm]
with Windows	KW-D1810	2.0	1000
	KW-D1816	2.0	1600
	KW-D2416T	2.8	1600
	KW-D2423T	2.8	2300
with Windows/ Spike	KH-D2416T	2.8	1600
	KH-D2423T	2.8	2300
Crocodile	KA-D2416T	2.8	1600
	KA-D2423T	2.8	2300
Crocodile with Spike	KB-D2416T	2.8	1600
	KB-D2423T	2.8	2300

## Single Use Retrieval Basket



Model Name	Compatible Endoscopes		Width of Basket [mm]
	Instrument Channel Diameter [mm]	Working Length [mm]	
GB-D1807L015	2.0	700	15
GB-D1810L020	2.0	1000	20
GB-D1819L020	2.0	1900	20
GB-D2415L030	2.8	1500	30
GB-D2419L020	2.8	1900	20
GB-D2423L035	2.8	2300	35

## Single Use Polypectomy Snare



Description	Model Name	Compatible Endoscopes		Width of Snare Wire [mm]
		Instrument Channel Diameter [mm]	Working Length [mm]	
Oval	DO-D2416-15	2.8	1650	15
	DO-D2423-15	2.8	2300	15
	DO-D2416-20	2.8	1650	20
	DO-D2423-20	2.8	2300	20
	DO-D2416-25	2.8	1650	25
	DO-D2423-25	2.8	2300	25
Hexagonal	DH-D2416-15	2.8	1650	15
	DH-D2423-15	2.8	2300	15
	DH-D2423-20	2.8	2300	20
	DH-D2416-25	2.8	1650	25
	DH-D2423-25	2.8	2300	25

## Single Use Injection Needle



Needle Size	Model Name	Compatible Endoscopes		Needle Length [mm]
		Instrument Channel Diameter [mm]	Working Length [mm]	
23G	NI-D1816-T2304	2.0	1650	4
	NI-D1816-T2305	2.0	1650	5
	NI-D1816-T2306	2.0	1650	6
	NI-D1816-T2308	2.0	1650	8
	NI-D2416-T2304	2.8	1650	4
	NI-D2416-T2306	2.8	1650	6
	NI-D2423-T2304	2.8	2300	4
	NI-D2423-T2305	2.8	2300	5
	NI-D2423-T2306	2.8	2300	6
	25G	NI-D1816-T2504	2.0	1650
NI-D1816-T2505		2.0	1650	5
NI-D2416-T2504		2.8	1650	4
NI-D2416-T2505		2.8	1650	5
NI-D2416-T2506		2.8	1650	6
NI-D2423-T2504		2.8	2300	4
NI-D2423-T2505		2.8	2300	5
NI-D2423-T2506		2.8	2300	6

## Single Use Sprayer



Model Name	Compatible Endoscopes	
	Instrument Channel Diameter [mm]	Working Length [mm]
TJ-D2418PB	2.8	1800
TJ-D2423PB	2.8	2300