

## J10 Series

# Video Naso-Pharyngo-Laryngoscopes

*Integrating High-Definition Imaging with Patient Safety and Comfort*

—  
**Stunning HD Image Quality**

—  
**Revolutionary Patient Comfort**

—  
**Advanced Hygiene Features**




# J10 Series

J10 Series Video Naso-Pharyngo-Laryngoscopes combine a vivid HD image with improved patient safety and comfort. When connected to the DEFINA EPK-3000 video processor, the VNL-J10 Series is an optimal solution for in-office assessment and therapeutic procedures.



## Stunning HD Image Quality

- Next-generation CCD imaging chip produces a sharp, clear, high-definition image supporting visualization of structure and detail
- Re-engineered clear distal end improves light distribution, resulting in a balanced and evenly illuminated image




| New clear distal end design for improved light distribution (on right)

## Revolutionary Patient Comfort

- At 2.4mm, the UltraSlim HD (VNL8-J10) features the smallest distal end diameter on a high-definition video naso-pharyngo-laryngoscope currently available on the market



| Offering a portfolio of naso-pharyngo-laryngoscopes for your varied needs



| The J10 Series Video Naso-Pharyngo-Laryngoscopes are compatible with i-SCAN™ image enhancement technology available in the DEFINA EPK-3000 Video Processor

## Advanced Hygiene Features

- The Therapeutic HD (VNL15-J10) utilizes innovative (single-use) suction control valves
- J10 Series is compatible with STERRAD NX and STERRAD 100NX<sup>1,2</sup>



| Single-use Suction Control Valve

## DEFINA Platform – J10 Series

	UltraSlim HD <b>VNL8-J10</b>	Slim HD <b>VNL11-J10</b>	Therapeutic HD <b>VNL15-J10</b>
Distal end diameter [mm]	2.4	3.5	4.8
Instrument channel diameter [mm]	–	–	2.0
Working length [mm]	300	300	300
Field of view [°]	80	80	80
Depth of field [mm]	3–50	3–50	3–50
Angulation up/down [°]	130/130	130/130	130/130

## PENTAX Medical

### Americas Headquarters

3 Paragon Drive  
Montvale, NJ 07645  
Phone +1 800 431 5880  
Fax +1 201 391 4189  
pentaxmedical.com/us

### Canadian Head Office

6715 Millcreek Drive, Unit 1  
Mississauga, ON L5N 5V2  
Phone +1 800 750 5558  
Fax +1 905 286 5571  
pentaxmedical.ca

**References:** 1. For a complete list of compatible reprocessing agents contact your PENTAX Medical representative 2. STERRAD is a registered trademark of Advanced Sterilization Products, division of Ethicon US, LLC, a Johnson & Johnson company.

Copyright 2018. All Rights Reserved. All company and product names and marks contained within are federally registered trademarks, trademarks or service marks of PENTAX Medical, Inc. Product specifications are subject to change without notice, and without any obligation on the part of the manufacturer.

**PENTAX**  
**MEDICAL**  
*Excellence in Focus*